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Instructions for use

Paste for fingers and hands

NOTE: The product is NOT a cream but a sticky paste that dries onto the skin. Therefore use in conjunction with a plaster or compress to prevent the paste from peeling off the affected area.

Important information before use: Make sure you have a plaster or compress available before beginning the application*. One plaster needed per application. Ensure that the plaster is large enough to cover the affected area.

1 Before application, wash and dry your hands and the affected area(s) on your hands.

2 Apply a generous amount of paste directly on to the affected area(s) using your fingers.

3 When treating cracks on fingers or hands ensure that the paste covers both the cracks and their outer edges. It is important to ensure that the paste fills the crack fully.

4 The paste's adhesive properties are activated with moisture. Therefore, after application of the paste, it is essential to moisten your finger and use it to gently smooth over the paste to seal it on your skin.

5 Cover the paste with a plaster or compress to prevent it from peeling off the affected area. Do not remove the plaster/compress before the next application.

6 Before re-applying, make sure to remove all of the previous paste application with plenty of warm water.

Warning: Do not use if allergic to any of the ingredients. Discontinue use if irritation or rash develops. Avoid contact with eyes. Keep out of reach of children. Store cool and dry. Up to two applications a day (morning and evening) is sufficient to achieve desirable results. You may increase or decrease the number of daily applications depending on the severity of your condition. You may continue using this product until the cracks and fissures no longer cause you discomfort, or disappear.

Manufacturer:

bonyf AG
Heiligkreuz 16
FL-9490 VADUZ
Liechtenstein - EUROPE
www.bonyf.com

INGREDIENTS: Calcium/Sodium PVM/MA Copolymer, Olea Europaea (Olive Fruit) Oil, Xanthan Gum, Alcohol, Cellulose Gum, Trihydroxystearin, Silica Dimethyl Silylate, Silica, Lecithin, Zinc Oxide, Paraffinum Liquidum, Citric Acid, Menthol, Sodium Lauryl Sulfate, Lactose Monohydrate, Petrolatum, Sodium Bicarbonate, Sodium Chloride, Potassium Caraoate, Citrus Medica Limonum Oil, Sodium Carbonate, Mentha, Arvensis Powder, Menthyl Lactate, Citral, Limonene.

 Made in Switzerland

*Plaster/Compress not supplied with the product.